

## Instability Repair

*Draw Tight Suture-Based  
Push-In Suture Anchors*

Parcus Medical has joined **Anika**



Surgical Technique Guide

# Instability Repair *Draw Tight Suture-Based Anchor*




step  
**1**

Position the drill guide on the prepared glenoid rim.




step  
**2**

Advance the 1.8mm Draw Tight Drill Bit until the positive stop contacts the drill guide.



step  
**3**

Insert the Draw Tight Anchor into the joint space. Once the tip of the anchor is aligned with the axis of the bone socket, advance the anchor into the socket using a mallet with small, short strikes. Insertion is complete when the handle of the driver contacts the end of the Drill Guide handle.



step  
**4**

Alternately pull the white Deployment Sutures to deploy the Draw Tight Anchor.


# Instability Repair *Draw Tight Suture-Based Anchor*



step  
**5**

Pull the sliding suture firmly to confirm deployment. Cut the Deployment Sutures.\*

*\* If additional suture is desired, the white deployment sutures may be used.*



step  
**6**

Shuttle the suture through the labrum.



step  
**7**


Reduce lesion and fix the tissue using preferred knots. Cut the suture tails above the knots.



step  
**8**

Repeat as necessary.

# Draw Tight *Suture-Based Anchor*



Draw Tight Suture-Based Anchors are recommended for use in both large and small-joint repairs. Made with UHMWPE and a natural PEEK-OPTIMA power tip.

## Features & Benefits

### Small insertion footprint & sub-cortical fixation

- Bone sparing

### Deployment Sutures

- Tactile confirmation the anchor has been successfully deployed

## DrawTight Suture-Based Anchors

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10565	1.8	-	Draw Tight, push-in, suture-based anchor, w/1, #2 suture (wht/blu)	UHMWPE	1

## Draw Tight Suture-Based Anchors w/Suture Tape

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
11148	1.8	-	Draw Tight, push-in, suture-based anchor, w/1 1.6mm suture tape (wht/blu)	UHMWPE	1

## Instrumentation for Draw Tight Suture-Based Anchors

Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10873U	-	Small Joint	1.8 Draw Tight, drill bit, w/positive stop	Unsterile	Reusable
10590U	-	Shoulder	1.8 Draw Tight, drill bit, w/positive stop	Unsterile	Reusable
10823	-	Small Joint	8-point slotted drill guide, w/beaver tail handle	Unsterile	Reusable
10446	-	Shoulder	8-point drill guide, w/trocar tip obturator	Unsterile	Reusable
10406	3.4	Shoulder	Blunt tip obturator <sup>1</sup>	Unsterile	Reusable
10771	3.3	Shoulder	Conical tip obturator, 1.8mm cannulation <sup>1</sup>	Unsterile	Reusable
10405	3.4	Shoulder	Trocar tip obturator <sup>1</sup>	Unsterile	Reusable

<sup>1</sup>Use w/10446 & 10330



Parcus Medical has joined Anika

**PARCUS**

**Parcus Medical, LLC.**

6423 Parkland Drive, Sarasota, FL  
1-941-755-7965 • www.parcusmedical.com

©2021 Parcus Medical, LLC. | QD 8095 Rev 1 02/2021  
Parcus is a registered trademark of Parcus Medical, LLC.

**ANIKA**